

## Homocystein S-Methyltransferase

Cat. No. EXWM-1697

Lot. No. (See product label)

### Einleitung

#### Beschreibung

Das Enzym verwendet S-Adenosyl-L-methionin als Methylspender weniger aktiv als S-Methyl-L-methionin.

#### Synonyme

S-Adenosylmethionin-Homocystein-Transmethylase; S-Methylmethionin-Homocystein-Transmethylase; Adenosylmethionin-Transmethylase; Methylmethionin:Homocystein-Methyltransferase; Adenosylmethionin:Homocystein-Methyltransferase; Homocystein-Methylase; Homocystein-Methyltransferase; Homocystein-Transmethylase; L-Homocystein S-Methyltransferase; S-Adenosyl-L-Methionin:L-Homocystein-Methyltransferase; S-Adenosylmethionin-Homocystein-Transmethylase; S-Adenosylmethionin:Homocystein-Methyltransferase

### Produktinformation

#### Form

Flüssigkeit oder lyophilisiertes Pulver

#### EC-Nummer

EC 2.1.1.10

#### CAS-Nummer

9012-40-2

#### Reaktion

S-Methyl-L-Methionin + L-Homocystein = 2 L-Methionin

#### Hinweise

Dieser Artikel erfordert eine Sonderanfertigung und die Lieferzeit beträgt zwischen 5-9 Wochen. Wir können nach Ihren Spezifikationen maßgeschneidert produzieren.

### Lager- und Versandinformation

#### Lagerung

Lagern Sie es kurzfristig bei +4 °C. Für die Langzeitlagerung lagern Sie es bei -20 °C~-80 °C.